

## Product Information

Catalog #    Description

### Premixed Sample Buffers

161-0737	2x Laemmli Sample Buffer, 30 ml
161-0747	4x Laemmli Sample Buffer, 10 ml
161-0710	2-Mercaptoethanol, 25 ml
161-0738	Native Sample Buffer, 30 ml
161-0739	Tricine Sample Buffer, 30 ml
161-0767	5x Nucleic Acid Sample Buffer, 10 ml
161-0768	TBE-Urea Sample Buffer, 30 ml
161-0763	IEF Sample Buffer, 30 ml
161-0764	Zymogram Sample Buffer, 30 ml

### Premixed Buffers

161-0732	10x Tris/Glycine/SDS, 1 L
161-0772	10x Tris/Glycine/SDS, 5 L
161-0734	10x Tris/Glycine, 1 L
161-0771	10x Tris/Glycine, 5 L

**Bio-Rad Laboratories, Inc.**

2000 Alfred Nobel Dr., Hercules, CA 94547 USA

510-741-1000

4006028 Rev J

---

## 2x Laemmli Sample Buffer

Catalog #161-0737

## 4x Laemmli Sample Buffer

Catalog #161-0747

**BIO-RAD**

## Introduction

Bio-Rad's Laemmli sample buffers are based on the method of Laemmli (1970). The use of Laemmli sample buffers ensures optimal band resolution when preparing proteins for SDS-PAGE with Tris-glycine-SDS running buffer.

## Formulations

<b>2x Laemmli sample buffer</b>	65.8 mM Tris-HCl, pH 6.8 26.3% (w/v) glycerol 2.1% SDS 0.01% bromophenol blue
<b>4x Laemmli sample buffer</b>	277.8 mM Tris-HCl, pH 6.8 44.4% (v/v) glycerol 4.4% LDS 0.02% bromophenol blue
<b>Storage</b>	Room temperature
<b>Shelf life</b>	2 years from date of manufacture

## Instructions for Use

### 1. Add Reducing Agent

To obtain a final 1x concentration of 355 mM 2-mercaptoethanol

**2x Laemmli sample buffer:** Add 50  $\mu$ l of 2-mercaptoethanol per 950  $\mu$ l.

**4x Laemmli sample buffer:** Add 100  $\mu$ l of 2-mercaptoethanol per 900  $\mu$ l.

Alternatively, add dithiothreitol (DTT or Cleland's reagent) to a final 1x concentration of 50 mM.

**Note:** For best results, do not store sample buffer with 2-mercaptoethanol.

### 2. Dilute Sample

**2x Laemmli sample buffer:** Dilute 1 part sample with 1 part 2x Laemmli sample buffer.

**4x Laemmli sample buffer:** Dilute 3 parts sample with 1 part 4x Laemmli sample buffer. More sample buffer can be added if necessary.

## Reference

Laemmli UK (1970). Cleavage of structural proteins during the assembly of the head of bacteriophage T4, Nature 227. 680–685.